


Ultra-sensitive Immune Profiling of Autoimmune Diseases with 200-Plex NULISAseq™ Inflammation Panel

Xiao-Jun Ma¹, Joanne Beer¹, Ishara Ariyapala¹, Qinyu Hao¹, Wei Feng¹, Dwight Kuo¹, Yiyuan Yin¹, Yuling Luo¹
¹Alamar Biosciences, Inc., Fremont, CA


Abstract

Cytokines and chemokines are critical components of the immune system and play important roles in autoimmune diseases. Comprehensive profiling of proteins in blood can provide deeper insights into the mechanisms underlying this highly complex and heterogeneous group of diseases. However, many of these proteins are present at very low concentrations in plasma, below the limit of detection of current immunoassays. We recently developed a novel automated multiplex immunoassay technology, NULISA™, capable of attomolar-level sensitivity, and a 200-plex inflammation-focused panel targeting 124 cytokines/chemokines and 80 other important inflammation and immune response-related proteins. We analyzed 21 plasma samples from patients with rheumatoid arthritis (n=5), Sjögren's syndrome (n=5), systemic lupus erythematosus (n=5), and ulcerative colitis (n=6), and 79 healthy donor samples using the 200-plex NULISAseq™ Inflammation Panel and the Olink Explore 384 Inflammation assay. Differential expression analysis using linear models identified 114 (56%) and 94 (26%) significant targets at a 5% false discovery rate in the NULISAseq and Olink Explore datasets, respectively. Among the 92 targets shared between these two platforms, 54 and 43 significant targets were identified by NULISAseq and Olink, respectively, with 36 in common. Many of the targets identified only by NULISAseq were low-abundance targets that were poorly detected by the Olink assay (IL4, IL5, IL20, IL17A, IL17F, IL33, and IL2RB) but have important roles in autoimmune diseases. In summary, with improved sensitivity and the most comprehensive inflammatory cytokine/chemokine panel, NULISA promises to be a powerful discovery tool for autoimmune disease research, which may lead to new diagnostic biomarkers and therapeutic targets.

NULISA Technology and Platform






Schematic of the NULISA workflow. 1) Immunocomplex formation; 2) first capture of immunocomplexes to dT beads; 3) bead washing to remove unbound antibodies and sample matrix components; 4) release of immunocomplexes into solution; 5) recapture of immunocomplexes onto streptavidin beads; 6) bead washing and DNA strand ligation to generate reporter DNA; 7a) detection and quantification of reporter DNA levels by qPCR; 7b) quantification of reporter DNA levels by NGS.






List of selected targets in 200-plex inflammation panel

AGER	AGRP	ANGPT1	C1QA	CCL13	CCL17	CCL20
CCL22	CCL23	CCL24	CCL25	CCL28	CCL3	CCL4
CCL7	CD200	CD200R1	CD276	CD40	CD70	CD83
CLEC4A	CSF1	CSF3	CST7	CXADR	CXCL1	CXCL10
CXCL12	CXCL14	CXCL3	CXCL8	CXCL9	EGF	EPO
FASLG	FGF19	FGF2	FLT3LG	GZMB	HGF	IFNG
IL10	IL10RB	IL11	IL12B	IL12RB1	IL13	IL15RA
IL16	IL17A	IL17C	IL17F	IL17RB	IL18	IL18R1
IL1B	IL1R2	IL1RN	IL20	IL2RB	IL32	IL33
IL3RA	IL4	IL4R	IL5	IL5RA	IL6	IL7
IRAK4	LAMP3	LGALS9	LTBR	MERTK	MMP10	NCR1
OSM	PTX3	TGFB1	TLR3	TNF	TNFRSF11A	TNFRSF13C
TNFRSF14	TNFRSF4	TNFSF10	TNFSF11	TNFSF12	TNFSF13	VEGFD

Results



200-plex NULISAseq™ Inflammation Panel performance characterization. **a.** The dynamic range for each target is indicated by the dark blue region. Values above the limit of detection (LoD) but below the lower limit of quantitation (LLOQ) are shown in lighter blue. Targets are ordered according to the geometric mean of the LLoQ and upper limit of quantitation (ULOQ). Dark blue lines represent generalized additive model (GAM) cubic regression spline fits to the ULoQ and LLoQ. The overall dynamic range for 200-plex NULISAseq spanned 9.6 log₁₀ values. **b.** Density plots of intraplate CV after internal control normalization (top) and interplate CV after internal control and intensity normalization (bottom). **c.** Cross reactivity. Two sets of 45 random antigen pools containing 4-5 targets each were analyzed with 200-plex NULISAseq. Each cell of the heatmap represents the percent of normalized read counts for that target (rows) occurring in that pool (columns). Targets were ordered according to pool membership such that the cells on the diagonal corresponded to the assigned pools. **d.** Detectability of 204 targets in 151 samples, including 79 from healthy controls and 72 from patients with various diseases. The y-axis represents NULISA Protein Quantification (NPQ) units minus the LoD for the respective target.



Conclusions

- NULISA achieves attomolar sensitivity by incorporating a unique background suppression mechanism to reduce assay background to a minimum
- NULISAseq enables high multiplexing and sample throughput while maintaining single attomolar sensitivity
- NULISA's ultrahigh sensitivity enabled detection low-abundance targets that were poorly detected by the Olink assay (e.g., IL4, IL5, IL20, IL17A, IL17F, IL33, and IL2RB) but have important roles in autoimmune diseases.
- NULISA provides a powerful new tool for both biomarker discovery and validation, which should facilitate seamless translational research.

Materials and Methods

Patient Samples

Human plasma samples were purchased from BioIVT (Westbury, NY).

Methods

21 plasma samples from patients with rheumatoid arthritis (n=5), Sjögren's syndrome (n=5), systemic lupus erythematosus (n=5), and ulcerative colitis (n=6), and 79 healthy donor samples were analyzed by using the 200-plex NULISAseq™ Inflammation Panel and the Olink Explore 384 Inflammation Panel. The NULISA assay was performed at Alamar Biosciences (Fremont, CA), and the Olink PEA assay was performed by Fulgent Genetics (Temple City, CA) and the High Throughput Biomarker Core of Vanderbilt University (Nashville, TN).

References:

Feng, W., Beer, J.C., Hao, Q. *et al.* NULISA: a proteomic liquid biopsy platform with attomolar sensitivity and high multiplexing. *Nat Commun* 14, 7238 (2023). <https://doi.org/10.1038/s41467-023-42834-x>



View Pub

